VETISEAL

(STERILE HEMOSTATIC AGENT DEVICE)

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◆ HEMOSTATA, INTERNAL, ABSORBABLE

VETIFOAM utilizes the absorbent properties of the raw material, bovine Gelatin, and the hemostatic properties of Kaolin, and acts to control bleeding from the tissue surface during surgical operations. The absorbent intracorporeal hemostatic product in contact with the bleeding area removes secretions and blood. It has a hemostatic effect by being absorbed and forming a physical barrier, and is gradually broken down and absorbed into the body before being discharged.

◆ PURPOSE OF USE

It is an absorbable material used to stop bleeding in body tissues during surgical operations.

◆ PREPARATIONS BEFORE USE

- This product is supplied sterilized by radiation (Gamma radiation, ISO11137-1,2,3), so re-sterilization is not permitted. The product packaging must be checked for damage and expiration date before use.
- Do not use if the packaging is damaged or the product is damaged.
- Before use, you must read the instructions carefully and observe all prohibitions, warnings, and precautions listed.
- Before use, the doctor must fully explain to the patient the contraindications, warnings, and precautions for use of this product and check the patient's condition.

◆ METHOD OF USE

Doctors must use it in accordance with its intended use and do not use it for purposes other than its intended use.

- Remove the final packaging aluminum pouch and the blister packaging inside it sequentially and take out the product.
- Insert the product into the area that requires hemostasis. Before insertion, adjust and use according to the affected area.
- Check the degree of hemostasis.

◆ MANAGEMENT METHOD AFTER USE

- This product is a disposable medical device, so reuse is prohibited after use and must be discarded.
- This product is transported in sterile packaging and must be undamaged upon arrival.

◆ STORAGE METHOD AND EXPIRY DATE

- It must be stored in a clean, dry place and at an appropriate temperature. (room temperature, 1 to 30 degrees)
- The expiration date of this product is 3 years from the date of manufacture.

◆ PRECAUTIONS FOR USE

- This product can be used by doctors who have been trained in its purpose and method of use, and must be used in accordance with the instructions provided by the manufacturer.
- The appropriate size and quantity must be selected by checking the degree of bleeding occurring in the patient's affected area.
- Be sure to check the expiration date of this product before using it. Do not use products whose expiration date has expired.

◆ CONTRAINDICATIONS

Taboos may be relative or absolute. The following may reduce the likelihood of a successful outcome:

- Patients with infection or contamination at the surgical site
 - Patients capable of producing an immune response
 - Patients allergic to gelatin protein
 - Patients who experience unexpected local and systemic allergic reactions
 - Persons judged by the surgeon to be unusable
 - Patients with bleeding or platelet disorders
 - Patients with unusual bleeding-related constitution or known coagulation disorders
 - Patients with autoimmune diseases
 - Patients with uncontrolled hypertension
- However, contraindications are not limited to this list.

◆ WARNING

- This product is disposable, so do not re-sterilize or reuse it. The manufacturer is not responsible for any damages or costs that may arise from reusing this product.
- Do not use this product if its packaging is damaged or opened before use.

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◆ REUSE

This product is a disposable sterilized veterinary medical device and should never be reused.

◆ HANDLING

In the event of side effects or accidents, it must be reported to the manufacturer.

◆ COMPLAINTS

Upon receiving any complaints or complaints regarding the quality, identity, durability, reliability, validity, and/or performance of the product, the user should notify manufacturer or their local representative. Furthermore, if the instrument is malfunctioning or suspected to be malfunctioning, it should be immediately notified to manufacturer or its local representative. If products of manufacturer have been improperly operated or have caused or contributed to the patient's serious injury or death, contact its local representative or manufacturer by phone, fax, or document as soon as possible.

Please fill out the following when receiving all complaints.

- Manufacturing number of components (Lot. Number)
- Product name (Part name)
- Part.No
- Name and address to be contacted
- Explain the cause of the complaint in detail

For additional information or complaints, please contact:

※Please check the attached label for product and manufacturer information.

Manufacturer : Mantiz Co., LTD.
Tel : +82 53-252-5833
Fax : +82 53-252-5834