



ABOVE THE STANDARD
MANTIZ



HEADQUARTER 423-15, Hyeoksin-daero, Dong-gu, Daegu, Korea 41071
TEL +82(0)53 252 5833 FAX +82(0)53 252 5834

SEOUL OFFICE #803, 17, Deogan-ro 104beon-gil, Gwangmyeong-si, Gyeonggi-do, Korea 14353
TEL +82(0)2 864 5833 FAX +82(0)2 864 5834

BUSAN OFFICE #912, 48, Centum jungang-ro, Haeundae-gu, Busan, Korea 48059
TEL +82(0)51 465 7117 FAX +82(0)51 465 7167

© 2024 MANTIZ Co., Ltd. All Rights Reserved
MZSC-H001-019 (Rev.000)

MANTIZ

Foam Series


Dressing, Absorbable, Wound, Biological Origin



Hemostatic Mechanism

Biocompatible materials ;


Development of optimal composition ratio of gelatin and kaolin for platelet aggregation



Mode of Action

<Hemostatic effect>

※ The role of each raw material in the hemostatic stage.



Features

POINT 01 | Complex Formulations

HEMOSTATIC EFFECT

'GELATIN+KAOLIN' Synergy effect

Bleeding control	Optimal composition ratio for platelet coagulation	Anti-inflammatory effect
------------------	--	--------------------------

POINT 02 | Absorbable

BIOCOMPATIBILITY

Use of Biocompatible Materials

Biodegradable bovine gelatin	Certified safe medical-grade kaolin
------------------------------	-------------------------------------

POINT 03 | Easy to use

USER CONVENIENCE

Simple and Easy

Product development reflecting user needs	Simple components and packaging	Various sizes and configurations
---	---------------------------------	----------------------------------



Product Information

Product	Classification	Ingredient	Specification	Expiry Date	Storage
			1.0 X 1.0 cm		
	Wound Dressing / Grade IV(MFDS)	Gelatin Kaolin	7.5 X 2.5 cm	1 years from Date of Manufacture	Room Temperature (1~25°C)
			3.75 X 3.75 cm 7.5 X 7.5cm		